


# LabSafe Series

## -25°C/Under-counter/ MDF-25U150

ANTECH undercounter lab freezer MDF-25U150 is ideal for small medical practices, pharmacies, and labs. Its height is designed to match laboratory table and optimized usage of limited space. Despite a compact freezer, the functions are complete for medical and laboratory applications, including data storage and USB port for download, high/low temperature alarm, door open alarm, power failure alarm, remote alarm contact, LED light etc. Therefore, MDF-25U150 meets the most rigorous requirements of healthcare and scientific research.

This unit uses environmentally friendly, natural refrigerants that are fully compliant with new EPA/SNAP/UL regulations.

MDF-25U150 is a stackable freezer (using Antech stacking kits) to install in a limited space. You can order with various combinations with flexibility for your specific demand.



▼ MDF-25U150

▼ With Stacking Kits



NATURAL REFRIGERANT



USB PORT



HIGH/LOW TEMPERATURE ALARM



DOOR AJAR A LARM



REVERSIBLE DOORS



STACKABLE DESIGN

MODEL	MDF-25U150
Temperature range (°C)	-10 ~ -25
Capacity (L)	120
External dimensions W×D×H (mm)	600×560×820
Shelves (pcs)	3
Data storage & USB port	•
Power supply	220V/50Hz; 220V/60Hz; 120V/60Hz