

SVF-100DR

For LyOrganic series freeze dryer

- Floating Scroll Technology
- 100% oil free
- Durable for organic solvent
- Maintenance free
- Compact size ideal for portable equipment
- Low running noise level



Technical Specification

Model	SVF-100DR
Medium	Oil-free dry gas
Pumping rate (L/min)	100
Ultimate vacuum (mbar)	≤0.05
Voltage (V/Hz)	220/50; 220/60; 120/60
Motor Power (W)	200
Inlet / Outlet Interface	KF 25 / G 3/8"
Ambient Temperature (°C)	-20~40
Noise [dB(A)]	≤55
Weight (kg)	8.5
Cooling system	Air Cooled
Dimensions (mm)	383*220*237